

MASSAGE CHAIR JP-4000
CYBER-RELAX AI

SPECIFICATIONS

Weight Approx. 96 kg
Power consumption 125 W

PLACEMENT & INSTALLATION

RECOMMENDED SPACE

To position the chair in your home, Ensure the doorway is at least 70 cm wide for the body and an additional 2 cm or more for carrying.


RECOMMENDED SPACE
72 cm WIDTH

INSTALLATION

Ensure the chair is installed on a flat surface and allocate sufficient space for reclining. Leave at least 30 cm for footrest extension and 50 cm for backrest reclining.



DIMENSIONS




* All dimensions are approximate.

IN 1954, THE WORLD'S FIRST
MASS-PRODUCED MASSAGE CHAIR
WAS DEVELOPED BY FUJIIRYOKI.

FUJIIRYOKI

FUJIIRYOKI HAS SOLD (SHIPPED)
MORE THAN 5.16 MILLION
MASSAGE CHAIRS.

THE HISTORY OF MASSAGE CHAIRS IS THE HISTORY OF FUJIIRYOKI



FUJIIRYOKI INDIA

A Division of Human Care International Pvt. Ltd.

Ph. : +91-9799979739 | Email : info@fujiiryoki.co.in

Website : www.fujiiryoki.co.in

CE CB RoHS ISO 13485:2016

Dealer

CYBER-RELAX

AI

MASSAGE CHAIR | JP-4000



FUJIIRYOKI

The design and specifications are subject to change without notice due to product improvement and other factors. Please note that there may be color differences between the print and the actual product. Photos and illustrations are for reference only. The contents of this catalogue are as of October 2024.

**EXPERIENCE
THE JAPANESE CRAFTSMANSHIP**

Massage therapists understand the unique body shapes and daily conditions of each individual, adjusting their techniques with meticulous care.

Fujiiryoki's Cyber-Relax AI aims to replicate these delicate human senses and movements, designed to deliver a massage that feels just right.

Enjoy the massage experience of a massage meister, effectively targeting and relieving tension points without being too strong or too weak.

Made in JAPAN



KEY HIGHLIGHTS

**CONDITIONING
MASSAGE**

Tailored to daily physical changes for total body care

5D-AI MECHA PLUS

With newly evolved massage techniques

40 UNIQUE COURSE MENUS

For a special experience every day

**13 VERSATILE STRETCH
PROGRAMS**

To stretch and relax muscle tensions of the whole body or specific body parts for a perfect tune-up

**8-INCH LARGE
TOUCH PANEL**

Provides intuitive control with Fujiiryoki's AI massage

**HIGH-FUNCTION
AIR SYSTEM**

With a carefully calculated airbag configuration and a fine strength adjustment function

DOUBLE HEATING

Gently warms you, enveloping you in comfort

TOTAL RELAX POSITION

For a highly relaxing reclining posture

5D-AI PLUS



5D-AI MECHA PLUS

WITH NEWLY EVOLVED MASSAGE TECHNIQUES



REAL-TIME AI MASSAGE ADJUSTMENT HIGH-PRECISION BRUSHLESS MOTOR


The AI provides delicate massage control and automatically adjusts the intensity in real-time according to the condition of your muscles, ensuring a stable massage that effectively targets and relieves tension points.



MASSAGES ARE PRECISELY TAILORED TO EACH INDIVIDUAL'S BODY


AI DOUBLE SENSING + (FUJIIRYOKI ORIGINAL)

Using real measurements, the "Shoulder Position Detection system" sets the height of the user's shoulders. Then, the "Backline Detection System" defines the actual spinal curve by measuring 8 different points, allowing the "5D-AI Mecha Plus" to personalize the pressure along the contours for the most comfortable massage.



EFFECTIVE NECK AND BACK MASSAGE WITH TWO MASSAGE BALLS

Our product's specialized two-ball massage system offers more precise and delicate movements compared to the four-ball systems. This allows for accurate targeting of specific points on the human body, ensuring a highly effective massage experience.



Adjusts the depth of the massage flexibly to match the contours of your body.

Moves while adjusting the width of the massage to conform to the shape of your shoulder blades and back.

Through the AI control system, It automatically reflects each user's body shape, Muscle status, and other information into the massage's intensity, techniques, and pace.

5D-AI MASSAGE MOVEMENT

HUMAN HAND-LIKE APPROACH



Massage balls gently contact the body to minimize the mechanical sensation.

THROUGH SHIATSU



When the massage balls protrude (Shiatsu), the movement gradually slows down, allowing it to carefully press into the body.

FIRM KNEADING




The AI automatically adapts to deliver a firm kneading massage on tight spots with consistent intensity.



98 MASSAGE TECHNIQUES INCLUDING 53 DIFFERENT KNEADING TECHNIQUES

The high-precision two massage ball mechanism delivers techniques tailored to the specific areas and conditions of your body, effectively relieving stiffness.




RHOMBOID

The two massage balls are designed to precisely fit the narrow area inside the shoulder blades, providing a kneading and relaxing massage. Ideal for relieving neck and shoulder fatigue from desk work.



3D WAVELET

This technique simultaneously performs "Tapping" and "Kneading," something that human hands cannot achieve. It provides a soothing relief for a tired lower back.



53 DIFFERENT KNEADING TECHNIQUES

KNEADING

KNEAD UP/DOWN
LOOP KNEAD UP/DOWN
SPIRAL KNEAD
PRECISION SPIRAL KNEAD
RELAX KNEAD UP/DOWN
PALMAR KNEAD UP/DOWN
3D KNEAD UP/DOWN
PUSH & KNEAD
RUBBING MASSAGE

TAPPING

TAPPING
WAVELET
IRREGULAR WAVES WAVELET
TAPPING PRECISION SPIRAL KNEAD
3D TAPPING
3D WAVELET

STRETCH

ROLLING
WAVE ROLLING
SHIATSU
SHIATSU STRETCH
STRETCH
STRETCH KNEAD UP
STRETCH TAPPING
STRETCH WAVELET
STRETCH SHIATSU
LIFT UP/PUSH DOWN

NECK/SHOULDER

NECK SHIATSU
NECK RELAX
NECK UPWARD LIFTING
NECK PALMAR KNEAD
NECK/SHOULDER EXTREME KNEAD
NECK/SHOULDER EXTREME TAPPING
SHOULDER TAPPING
UPPER SHOULDER PRESS
COLLAR DOWNWARD PRESSING
TOP SHOULDER LINE SHIATSU


BACK/WAIST

SHOULDER BLADE LINE
RHOMBOID
SHOULDER BLADE OPEN
UPPER BACK DEEP SHIATSU
WAIST EXTREME KNEAD
WAIST EXTREME TAPPING

PARASPINAL
LUMBAR DEEP SHIATSU
WAIST SHIATSU
WAIST PALMAR
BUTTOCK


13 DIVERSE STRETCHING PROGRAMS

With a full-flat position, you can thoroughly stretch your hip joints. It helps to extend your contracted body while promoting blood circulation through massage.



**FULL BODY STRETCH:
6 COURSES**

**TARGETED STRETCH:
7 COURSES**



LARGE 8-INCH TOUCH PANEL

Unlock the full potential of Fujiiryoki's AI technology with our intuitive large touch panel, allowing you to effortlessly control your massage experience with just a single touch.




FOOTREST PERFECT FIT

The double sliding mechanism (electric and spring) allows to accommodate various users heights and leg types.

THOROUGH FULL-BODY MASSAGE WITH HIGH-FUNCTION AIR SYSTEM

The airbags envelop your entire body, refreshing your accumulated fatigue. Each of 8 air massage regions can be adjusted in 7 strength levels. You can experience rhythmic movements in pulse mode that promote blood circulation, as well as an air massage with varying speeds that mimics a real hand massage.




AIRBAGS
8 BODY REGIONS
7 LEVELS INTENSITY
ADJUSTMENT

HIGH-SPEED INFLATION AIRBAG SYSTEM

The built-in high-flow and high-pressure air pump allows the high-capacity airbags to operate powerfully, enhancing the massage effect.

ADJUSTABLE ARMREST

The armrests can be slid to fit the length of your arms, Ensuring a snug fit even when reclining.



FOOT SOLE SHIATSU

After securely gripping the ankles and insteps with air bags, Dedicated rollers provide focused massage to the entire surface of the feet.

**Shiatsu* is here used to describe the roller movements

40 UNIQUE COURSE MENUS


Experience 40 expertly designed courses, each crafted under the supervision of a massage master to provide an ideal massage tailored to your current needs.

FULL BODY
9 COURSES

TARGETED AREAS
15 COURSES


STRETCH
13 COURSES

RECOMMENDED
2 COURSES




EASY SELECTION

Quickly choose from two newly recommended courses with simple controls.



COURSE RECOMMENDATION

Get the best course recommendations based on your fatigue level.




PURPOSE/SITUATION-BASED SELECTION

Smoothly select courses based on your specific goals or situations.


USER CUSTOMIZE MODE

Register your preferred massage settings for a customized experience, and your exclusive massage will begin immediately. You can finely adjust the Shoulder position, Spinal line, Air pressure, and Reclining angle, with options to save settings for up to 2 users.



BODY SHAPE MEASUREMENT AI DOUBLE SENSING+

Thanks to the "shoulder position detection system", which identifies the shoulder height, and the "s-line detection system", which detects spinal alignment, the AI can automatically optimize the massage points.



REGISTRATION OF MASSAGE INTENSITY FROM SHOULDERS TO HIPS CUSTOMIZABLE SPINAL LINE

Based on the detected back line, the consistent massage intensity is applied from the shoulders to the buttocks. In addition, the kneading intensity of each of the 8 spinal points can be adjusted according to user preferences.



REGISTRATION OF MOST RELAXING POSITION CUSTOMIZABLE RECLINING ANGLE

The reclining angle can be freely adjusted and registered.



ADJUSTING THE STRENGTH OF THE AIR PRESSURE THROUGHOUT THE BODY CUSTOMIZABLE AIR PRESSURE

You can customize the strength of each of the 8 different air massage sections by experiencing the real-time effects and registering your preferences for future use.



BLUETOOTH SPEAKER

Surround yourself with your favorite sounds and create an even more relaxing atmosphere.



MINDFULNESS COURSE [AUTOMATIC]

Based on the concept of "meditation", a 30-minute massage performed through airbags only. The slow air massage carefully and gently soothes the entire body.

DOUBLE HEATING

Two heating systems are built in to create a more comfortable massage time.



QUICK- HEAT MULTI-HEATERS

The pads with built-in warming function can be placed behind the back to warm up around the waist, or can be unzipped and placed on the front to warm around the abdomen.



FOOT SOLE HEATING

Gently warming up the feet, which tend to get cold easily.

1. The heating function can be used independently.
2. Set temperature: about 38-45C° [at ambient temperature 20-25C°]

* Set temperature may feel different from perceived temperature.



TOTAL RELAXING POSITION

The relaxation level is enhanced by the comfortable lying position, which is also excellent for optimizing the stretching effects.

A remote control stand that allows for easy touch operation, even while lying down.



ONYX(BK)



IVORY(CS)



CHOCOLATE(BG)