

Massage Chair MR8000  
CYBER-RELAX PRO


Specifications

Weight Approx. 148kg  
Power consumption 210W

Placement & Installation

Recommended Space


To position the chair in your home, ensure the doorway is at least 76 cm wide for the body and an additional 2 cm or more for carrying.




Recommended Space  
78 cm Width

Installation

Ensure the chair is installed on a flat surface and allocate sufficient space for reclining. Leave at least 22 cm for footrest extension and 10 cm for backrest reclining.



Dimensions



\*All dimensions are approximate.

IN 1954, THE WORLD'S FIRST  
MASS-PRODUCED MASSAGE CHAIR  
WAS DEVELOPED BY FUJIIRYOKI.

FUJIIRYOKI

FUJIIRYOKI HAS SOLD (SHIPPED)  
MORE THAN 5.16 MILLION  
MASSAGE CHAIRS.

THE HISTORY OF MASSAGE CHAIRS IS THE HISTORY OF FUJIIRYOKI

1954



1954 The first  
Fuji automatic  
massage chair.

1967 Fuji C-3  
Specializes in  
tapping massage.



1967 Fuji FM-52  
Massage can be  
changed with a  
single switch.



1987 Fuji KIT-10DX  
Specializes in  
rolling massage.



1987 Fuji MC-133  
World's first air-based  
massage chair.



1995 Fuji Cyber-Relax  
AS-2100  
Premium deep-tissue  
massage.



2001 Fuji Cyber-Relax  
AS-2100  
Premium deep-tissue  
massage.



2020 Fuji Cyber Relax AS-001  
A smart massage chair.



Massage Chair MR8000  
CYBER-RELAX PRO  
High-function next-generation  
massage chair.

\*As of December 31, 2023  
from 1965/ Researched by  
FUJI MEDICAL INSTRUMENTS MFG. CO., LTD.

2024

FUJIIRYOKI INDIA

A Division of Human Care International Pvt. Ltd.

Ph. : +91-9799979739 | Email : info@fujiiiryoki.co.in

Website : www.fujiiiryoki.co.in

CE CB RoHS ISO 13485:2016

Dealer

Blank box for dealer information.

CYBER-RELAX

PRO

MASSAGE CHAIR | MR8000



FUJIIRYOKI

THE DESIGN AND SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE DUE TO PRODUCT IMPROVEMENT AND OTHER FACTORS.  
PLEASE NOTE THAT THERE MAY BE COLOR DIFFERENCES BETWEEN THE PRINT AND THE ACTUAL PRODUCT.  
PHOTOS AND ILLUSTRATIONS ARE FOR REFERENCE ONLY. THE CONTENTS OF THIS CATALOGUE ARE AS OF APRIL 2025.



# RIDE TO THE FUTURE

## Next-Generation High-Function Massage Chair

Elevate your tomorrow.

A serene oasis to recharge your mind and body amidst the bustle of your day.

CYBER-RELAX PRO delivers effective, full-body care in no time, helping you conquer fatigue and embrace comfort.

Discover our massage chair, crafted to ensure you greet tomorrow in peak condition.

FLEX SL TRACK

DUAL MESSAGE  
ROBOTS

19 DISTINCT  
MESSAGE TECHNIQUES

HIGH PERFORMANCE  
AIR MESSAGE SYSTEM

CALF ROLLER /  
TRIPLE FOOT ROLLER

LAY FLAT STRETCH

LUXURIOUS  
FULL-BODY HEATER

VARIETY OF  
20 AUTO COURSES

Designed by FUJIIRYOKI in JAPAN



## FLEX SL TRACK

Equipped with the next-generation global standard Flex SL frame, this variable frame allows for both the previously unattainable Zero Gravity posture and the Lay Flat posture, offering the ultimate in relaxation.

## THE 130CM-LONG RAIL COVERS FROM NECK TO THIGHS

The massage robot moves smoothly along the variable sections while the chair reclines at any angle, providing a gentle yet thorough massage for the entire body.

### UPRIGHT



It can be installed in limited spaces, as it requires only about 10 cm of rear clearance to fully recline.

### ZERO-GRAVITY



The body is tilted downward, with the legs slightly raised to support the spine and alleviate pressure on the discs in the back, allowing you to relax during the massage.

### LAY-FLAT




This position is ideal for stretching the body in a flat position while being

## DUAL MASSAGE ROBOTS


Super-efficient massage system that uses two dedicated mechanisms to massage the entire body at once.

With two mechanisms moving independently along the rails, not only does it provide total relaxation care for your upper and lower body parts, but it can also provide a customized massage targeting specific muscle groups as intended, leaving you feeling massaged by multiple masseurs.




### 4D QUAD-BALL MECHANISM

The high-function mechanism with a maximum protrusion of 7 cm and an actual vertical travel area of 65 cm smoothly relaxes the neck and entire upper body.



### 2D QUAD-BALL MECHANISM

The full-fledged mechanism, which operates over a wide area centered on the seat surface with an actual vertical travel area of 62 cm, dynamically relaxes the lower body.



### FULL BODY MASSAGE



Each mechanism is divided into upper and lower sections to efficiently cover the entire body.

### UPPER BODY MASSAGE




Both mechanisms are fine-tuned around the shoulders/back with a narrower range on the upper side.

### LOWER BODY MASSAGE



Two mechanisms focus on the lower side of the body, intensively and dynamically massaging the lower back and hips.





## MASSAGE TECHNIQUES

19 distinct massage techniques mimic the natural movements of human hands.

Aiming to replicate complex hand movements by combining techniques such as Kneading, Tapping, and Shiatsu. Numerous techniques are incorporated to approach muscles and stiffness in an optimal manner for effective massage and relaxation.

KNEAD UP	KNEAD DOWN	4D KNEAD UP	4D KNEAD DOWN
TAPPING	4D TAPPING	WAVELET	4D WAVELET
ROLLING	SHIATSU	STRETCH SHIATSU	DEEP SHIATSU
SHOULDER GRASP	3D KNEAD UP	3D KNEAD DOWN	3D TAPPING
3D WAVELET	PALM KNEAD UP	PALM KNEAD DOWN	

### MIX

Conventionally, only 1 massage technique was selected for manual massage. However, by mixing up to 4 massage techniques, even a manual massage can transcend simple repetition, allowing more than 1,000 combinations and realizing a more diverse customized experience.

\*"Shiatsu" is a massage action in which vertical pressure is applied to muscles with a knead ball.

## HIGH PERFORMANCE AIR MASSAGE SYSTEM

Featuring 86 strategically placed air cells, this whole-body air massage system promotes highly efficient blood flow.

Airbags are meticulously placed throughout the entire body, enveloping it in comfort. A single arm utilizes 20 individual airbags, mimicking the precise movements of human fingers. This detailed approach stimulates blood flow and provides a comprehensive


### NECK-SPECIFIC AIRBAG MASSAGE & HEATER

Experience a revolutionary neck massage that goes beyond traditional kneading with massage balls.




## CALF ROLLER / TRIPLE FOOT ROLLER

Airbags and rollers are luxuriously used. The roller massages the foot while the airbag provides firm pressure, promoting blood circulation and relieving fatigue.



### CALF

Each leg is wrapped with five airbags: two on each side and one in the back on each side of the calf. A special calf roller with protrusions relaxes the muscles in the stiffened



### FOOT

A large air bag holds the ankle and heel firmly in place, while triple foot rollers apply acupressure and massage the ankle





## VARIETY OF 20 AUTO COURSES

Original courses combining the massage techniques, airbags, and more, tailored to your physical condition and mood.

### REJUVENATE

- Swedish
- Revive
- Restore
- Shiatsu
- Energize

### STRETCH

- Thai
- Dynamic
- Static
- Deep
- Relax

### FOCUS

- Neck
- Upper Back
- Lower Back
- Legs
- Glutes

### RELAX

- Unwind
- Calm
- Dream
- Comfort
- Meditation

### DIY COURSE

Create your own massage course. As you gain more experience with the massage chair, your preferences will naturally evolve. This new feature allows you to create your own massage programs to your evolving needs and personal tastes.

### LAY FLAT STRETCH

Full-body stretch movement utilizing the flexible track to stretch hip joints, thighs, and legs.

## LUXURIOUS FULL-BODY HEATER

Heater for Neck, Hands, Back/Stomach, and Calves.  
The whole body is wrapped in warmth with the rich heating system.



It is also possible to warm the stomach area by placing the heater pad from the back pad onto the front of the body.



## RICH SOUND SYSTEM

Immerse yourself in a healing space enveloped by high-quality speakers. Experience rich sound with dedicated high- and low-frequency speakers (woofer) for a truly immersive audio experience.

Connect wirelessly via Bluetooth and enjoy your favorite music or relax to the seven built-in relaxing sounds.





8-INCH TOUCH PANEL & STYLISH QUICK KEYS



8-INCH TOUCH PANEL

The user-friendly interface, featuring a large screen and intuitive menus, makes navigation a breeze. Customize your massage experience with detailed settings that are



STYLISH QUICK KEYS


Stress-free, quick access to basic menus such as on/off, reclining and favorite course etc.



GRAPHITE (DG)



MIDNIGHT (BK)



ADMIRAL (DB)



WHITE (WH)



CHAMPAGNE (LB)

WELCOME LIGHTING SYSTEM

MOTION SENSOR LIGHT

Sensors mounted on the massage chair react when a person approaches, and the side panel and head light glow as if the chair is welcoming you.



CONTACTLESS CHARGING PORT



Conveniently charge your phone wirelessly on the pad or with a USB cable while having a massage.

HEALTH MONITORING BY FINGER SENSING



Sensing biometric information from the fingertip and analyzing fatigue levels. This system then recommends the most suitable massage